



LDPE TRANSFER PIPETS – ACTUAL SIZE




GENERAL PURPOSE




Item #	Technical Data	Case Qty
AKM-3800-0051	Volume: 4ml Drops/ml: 25 Bulb Draw: 2ml	Total Length: 130mm (5.1") Stem Diameter: 6.6mm (0.26") Tip Diameter: 2.5mm (0.1")



Item #	Technical Data	Case Qty
AKM-3800-0037	Volume: 7ml Drops/ml: 20 Bulb Draw: 3.2ml	Total Length: 155mm (6.1") Stem Diameter: 7.8mm (0.31") Tip Diameter: 3.1mm (0.12")



Item #	Technical Data	Case Qty
AKM-3800-0052	Volume: 8ml Drops/ml: 20 Bulb Draw: 4.6ml	Total Length: 152mm (6.0") Stem Diameter: 6.2mm (0.24") Tip Diameter: 3mm (0.12")




Item #	Technical Data	Case Qty
238718	Volume: 8.5ml Drops/ml: 18 Bulb Draw: 4.9ml	Total Length: 137mm (5.4") Stem Diameter: 6.1mm (0.24") Tip Diameter: 4mm (0.16")


LDPE TRANSFER PIPETS – ACTUAL SIZE




GRADUATED




Item #	Technical Data	Case Qty	
238456	Volume: 1.5ml Drops/ml: 21 Bulb Draw: 0.9ml Total Length: 116mm (4.6")	Stem Diameter: 6.1mm (0.24") Tip Diameter: 4mm (0.16") Graduated to: 0.3ml Graduation Increments: 0.1ml	5,000



Item #	Technical Data	Case Qty	
AKM-3800-0033	Volume: 3ml Drops/ml: 25 Bulb Draw: 2.2ml Total Length: 140mm (5.5")	Stem Diameter: 4.4mm (0.17") Tip Diameter: 2.4mm (0.09") Graduated to: 1ml Graduation Increments: 0.25ml	5,000




Item #	Technical Data	Case Qty	
238692	Volume: 5ml Drops/ml: 23 Bulb Draw: 1.8ml Total Length: 155mm (6.1")	Stem Diameter: 6.5mm (0.25") Tip Diameter: 2.7mm (0.11") Graduated to: 2ml Graduation Increments: 0.5ml	5,000




Item #	Technical Data	Case Qty	
AKM-3800-0022	Volume: 5ml Drops/ml: 24 Bulb Draw: 3.4ml Total Length: 145mm (5.7")	Stem Diameter: 4.8mm (0.19") Tip Diameter: 2.7mm (0.11") Graduated to: 1ml Graduation Increments: 0.25ml	5,000

GRADUATED con't




Item #	Technical Data	Case Qty
AKM-3800-0023	Volume: 7ml Drops/ml: 20 Bulb Draw: 3.2ml Total Length: 155mm (6.1") Stem Diameter: 7.9mm (0.31") Tip Diameter: 3mm (0.12") Graduated to: 3ml Graduation Increments: 0.5ml	5,000



Item #	Technical Data	Case Qty
238716	Volume: 7ml Drops/ml: 21 Bulb Draw: 3.5ml Total Length: 148mm (5.8") Stem Diameter: 7.5mm (0.30") Tip Diameter: 3.2mm (0.13") Graduated to: 3ml Graduation Increments: 0.5ml	5,000

FINE TIP



Item #	Technical Data	Case Qty
238455	Volume: 1.3ml Drops/ml: 45 Bulb Draw: 1ml Total Length: 51mm (2.0") Stem Diameter: 1.15mm (0.05") Tip Diameter: 1.15mm (0.05")	4,000




Item #	Technical Data	Case Qty
AKM-3800-0034	Volume: 1.5ml Drops/ml: 50 Bulb Draw: 1ml Total Length: 104mm (4.1") Stem Diameter: 4.5mm (0.18") Tip Diameter: 0.9mm (0.04")	3,200


LDPE TRANSFER PIPETS – ACTUAL SIZE




FINE TIP con't




Item #	Technical Data	Case Qty
AKM-3800-0035	Volume: 5ml Drops/ml: 50 Bulb Draw: 3.2ml	Total Length: 153mm (4.1") Stem Diameter: 4.6mm (0.18") Tip Diameter: .85mm (0.03") 2,000



Item #	Technical Data	Case Qty
238666	Volume: 5.8ml Drops/ml: 50 Bulb Draw: 3.4ml	Total Length: 157mm (6.2") Stem Diameter: 4.6mm (0.18") Tip Diameter: .85mm (0.03") 5,000



Item #	Technical Data	Case Qty
238665	Volume: 5.8ml Drops/ml: 60 Bulb Draw: 3.4ml	Total Length: 147mm (5.8") Stem Diameter: 4.6mm (0.18") Tip Diameter: .85mm (0.03") 5,000




Item #	Technical Data	Case Qty
AKM-3800-0030	Volume: 8.7ml Drops/ml: 45 Bulb Draw: 4.8ml	Total Length: 147mm (5.8") Stem Diameter: 6.1mm (0.24") Tip Diameter: 1.4mm (0.06") 4,000


LDPE TRANSFER PIPETS – ACTUAL SIZE




NARROW STEM




Item #	Technical Data	Case Qty
238454	Volume: 1.2ml Drops/ml: 25 Bulb Draw:0.9ml	Total Length: 65mm (2.56") Stem Diameter: 2.5mm (0.1") Tip Diameter: 2.5mm (0.1") 5,000



Item #	Technical Data	Case Qty
238663	Volume: 4ml Drops/ml: 24 Bulb Draw:3ml	Total Length: 85mm (3.4") Stem Diameter: 2.5mm (0.1") Tip Diameter: 2.5mm (0.1") 5,000




Item #	Technical Data	Case Qty
AKM-3800-0032	Volume: 4.5ml Drops/ml: 25 Bulb Draw:3.2ml	Total Length: 155mm (6.1") Stem Diameter: 2.5mm (0.1") Tip Diameter: 2.5mm (0.1") 5,000




Item #	Technical Data	Case Qty
238703	Volume: 6.5ml Drops/ml: 28 Bulb Draw:3.2ml	Total Length: 155mm (6.1") Stem Diameter: 7.8mm (0.31") Tip Diameter: 2.1mm (0.08") 5,000

NARROW STEM con't




Item #	Technical Data		Case Qty
AKM-3800-0053	Volume: 15ml Drops/ml: 19 Bulb Draw: 8ml	Total Length: 155mm (6.1") Stem Diameter: 3mm (0.12") Tip Diameter: 3mm (0.12")	2,500

PADDLE PIPETS



Item #	Technical Data		Case Qty
238452	Volume: 0.8ml Drops/ml: 25 Bulb Draw: 0.2ml	Total Length: 125mm (4.9") Stem Diameter: 3.6mm (0.14") Tip Diameter: 3.6mm (0.14")	5,000



Item #	Technical Data		Case Qty
238552	Volume: 1ml Drops/ml: 38 Bulb Draw: 0.2ml	Total Length: 130mm (5.1") Stem Diameter: 3.8mm (0.15") Tip Diameter: 1.85mm (0.07")	5,000